

Guidance to Emergency Medical Services (EMS), Police Forces, Fire Departments, Skilled Nursing Facilities, Long Term Care Facilities, Residential Treatment Facilities, Mental Health Facilities, Dental Offices, Elder Care Facilities, Correctional Facilities, Visiting Nurses, and Home Health Aides on

N-95 Respirator Disinfection in Your Region

The State Emergency Operations Center (SEOC) has provided multiple hospitals with Steris V-Pro low temperature hydrogen peroxide sterilizing units. These units are typically used for sterilizing other medical devices but were recently granted an Emergency Use Authorization by the U.S. Food and Drug Administration to sterilize certain N-95 respirators.

Five hospitals in Vermont have received these units and will provide N-95 sterilization services for you during the COVID-19 Response. The five hospitals with points of contact are:

- University of Vermont Medical Center
 - Drew Robinson; Andrew.Robinson@uvmhealth.org; 802-847-4258
- Southwestern Vermont Medical Center
 - Joy Bull; Joy.Bull@svhealthcare.org; 802-442-6321 or Tanya Crowder; Tanya.Crowder@svhealthcare.org; 802-442-6361
- Brattleboro Memorial Hospital
 - Mike Geissler; mgeissler@bmvht.org; 802-257-8288
- Springfield Hospital
 - Jim Smith; JimSmith@springfieldmed.org; 802-885-7651
- Northeastern Vermont Regional Hospital
 - Pete Tomczyk; p.tomczyk@nvrh.org; 802-748-7348

While other hospitals may have similar sterilizer capabilities, they are not obligated to sterilize outside agency N-95 respirators like this, though they may choose to do so independently.

You should reach out to your region's hospital from the above list soon to discuss. Email may be better to reduce the impact. ***They will have requirements that must be met as you package your respirators for sterilization for you.*** We have asked them to diligently work with you during the COVID-19 Response for N-95 respirator sterilization. After COVID-19 Response ends, this obligation to sterilize for you will end, too.

Respirator Sterilization Processing

Respirators for sterilization must be very clean – **no writing, no lipstick or other makeup**. Tested and acceptable respirator types to date include (more will follow):

- 3M 8000, 8210, 1860, 1860S, 1870P
- Moldex 1510Z, 1511, 1512, 1513, 1517

Discard soiled or damaged respirators in accordance with normal procedures – if in doubt, throw it out; **keep only used respirators that are in good condition**. The sterilizing facility may not return all of your

respirators if they deem they are too soiled for sterilization or damaged for reuse. N-95 sterilization may occur a finite number of times.

Manage used, unsterilized respirators for disinfection as biohazardous materials. They must be packaged in accordance with OSHA Bloodborne Pathogen Standard and transported in accordance with applicable US DOT regulations. It is hoped that local licensed biosafety transportation in accordance with DOT regulations can be arranged by the you and the hospital, but the State can help with this during COVID-19 Response by calling the SEOC call takers at 800-347-0488.

Put used respirators together into a biohazard labeled plastic bag, line a box with a large plastic bag and put the bagged respirators into the bag lining the biohazard labeled cardboard box. The box should be labeled with the facility name and point of contact.

Please recall that N-95 respirator usage has three regulatory requirements for each user before use: medical evaluation for safety purposes, respirator training on safe usage, and fit-testing to ensure the size and type respirator affords protection for the individual. There are resources for each available locally and through the SEOC during the COVID-19 Response.

Should you have any questions that cannot be answered working locally, please contact the SEOC call takers at 800-347-0488.

Thank you for all you have done and are doing for the COVID-19 Response.

Vermont State Emergency Operations Center