

Data Logger Tips

- Download information **weekly** because the data logger has limited storage memory and when memory is full, the data logger will stop recording.
- You are still required to check and record temperatures **Max and Min once a day** and document time and initials for each reading.
- Current temperature is always displayed on the data logger.
- You should check and review Min/Max from previous day in the morning by pressing the **“Review”** button. Refer to Quick reference guide for more detailed information.
- Any time you press the “Review” button the data logger is paused for 30 minutes.
- Always pause the data logger before you disconnect it. If it takes you longer than 30 minutes to plug it in the cradle, the data logger will display -199°C or -40°C. By downloading the information from data logger -199°C or -40°C will be cleared.
- If the data logger is restarted and the cable is not all the way plugged in, the data logger will display inaccurate reading (may display -15°C on the refrigerator or colder). Download the data logger, plug in the cable all the way and restart it.
- Every time you finish a download, you must start the data logger by pressing and holding **“Start”** button until the word, “Starting,” stops flashing.
- Call the Immunization Program when:
 - * Current temperature is out of range
 - * Min/Max is out of range
 - * “Alarm” is displayed on data logger
- When **“Alarm”** displays on the data logger, download the information to clear the data logger.
- Use Immunization Program supplied temperature logs
- Technical support for data loggers is provided by Control Solutions Inc. by calling 1-888-311-0636
- For more information on downloading or reviewing refer to Quick Reference Guides provided during the data logger installation