

# How to Conduct the Intrivo On/Go COVID-19 Antigen Self-Test

**STEP 1.** Obtain the test from the testing site.

**STEP 2.** Review all instructions associated with the test.

**STEP 3.** Conduct the test off-site and wait 10-15 minutes for your results.

**STEP 4.** Report results at [healthvermont.gov/reportresults](https://healthvermont.gov/reportresults) and complete the short form.

It's important to report your test results.  
Your response is confidential.



**STEP 5.** If proof of your test results is required, download the On/Go App.



For translated testing instructions and instructional videos, please visit [www.healthvermont.gov/covid-19/testing](https://www.healthvermont.gov/covid-19/testing)

# FAQ

## **What is an Intrivo On/Go COVID-19 Antigen Self-Test?**

The Intrivo On/Go COVID-19 Antigen Self-Test is a lateral flow immunoassay intended for the qualitative detection of nucleocapsid protein antigen from SARS-CoV-2 from individuals with or without symptoms or other epidemiological reasons to suspect COVID-19 infection when tested twice over three days with at least 36 hours between tests. This test is authorized for non-prescription home use with self-collected direct anterior nasal (nares) swab samples from individuals aged 14 years or older, or adult-collected anterior nasal swab samples from individuals aged two years or older.

**How does a COVID-19 Antigen Self-Test compare to a molecular test?** Antigen tests are very specific for the virus, but are not as sensitive as molecular tests. This means that a positive result tends to be accurate, but a negative result does not rule out infection.

**When does my test expire?** Currently, the test expires on the last date of the month printed on your box (i.e. January 2022 means January 31, 2022). However, with medical products under EUA, such as On/Go, the FDA, after close studies, grants rolling expiration extensions. Visit this website for further details and see full list of extensions [www.letsongo.com/product-updates](http://www.letsongo.com/product-updates)

**Why are there two tests in the kit?** Per On/Go, the test is enabled for over-the-counter/self-use when people test themselves twice over three days with at least 24 hours (and no more than 48 hours) between tests. The two tests included allow individuals to fulfill the recommended use.

**What if my test results are positive?** If you test positive, be sure to report your results via the link on the other side of this card. For additional guidance, visit [healthvermont.gov/covid19positive](http://healthvermont.gov/covid19positive)



To learn more about the Intrivo On/GO COVID-19 Antigen Self-Test, visit [www.letsongo.com/faq](http://www.letsongo.com/faq)